
# Cell culture solutions

High quality ingredients for research, manufacturing and clinical applications

cQrex™







## Our offering: Boosters and beyond









### **BOOSTERS**

cGMP-quality peptides and keto acids solving key issues in mammalian cell culture

### **AMINO ACIDS**

More than two decades of excellence in the supply of high-quality, animalorigin free ingredients for consistent media performance

### **CUSTOMIZATION**

Media ingredients and formulations tailored to your specific requirements

#### **SERVICES**

Screenings and booster development support with strong local market expertise

### The cQrex<sup>TM</sup> effect: How our boosters work

Tyrosine and glutamine are examples of amino acids that are hardly soluble or chemically reactive and unstable. This limits process performance, especially if highly concentrated feed or perfusion media are required to intensify the process.

cQrex<sup>TM</sup> dipeptides and performance boosters from Evonik provide an elegant solution. Chemical coupling to carrier amino acids yields dipeptides that are stable and highly soluble at neutral pH. These natural building blocks are efficiently utilized by a variety of animal and human cells.

### **PRODUCT PROPERTIES**

- · High solubility in neutral pH
- · High stability
- · High bioavailability

### **PROCESS BENEFITS**

- Increased productivity
- No need for high-pH feeds, reduced production risks
- Process simplification
- Lower cost per unit of product

### **DESIGNED FOR CELL CULTURE**

## Curious? Order your sample kit!

When performance is of the essence, it helps to have the right tools available. The cQrex<sup>TM</sup> standard booster portfolio consists of peptides, keto acids, special amino acids and derivatives. They are designed to stabilize

components, to prolong longevity, solubilize critical constituents, increase productivity, and help manage the concentration of critical metabolites. All boosters are readily available at commercial scale.

PRODUCT NAME	BENEFITS	SOLUBILITY IN WATER (20°C, <sub>P</sub> H 7)	QUANTITY IN SAMPLE KIT
<b>cQrex™ GY</b> Glycyl-L-tyrosine	Tyrosine replacement when high concentrations at neutral pH are needed	34g/L	10 g
<b>cQrex™ GQ</b> Glycyl-L-glutamine	Alternative stable replacement of free glutamine	154g/L	10 g
<b>cQrex™ AQ</b> L-alanyl-L-glutamine	The classical stabilized form of glutamine	720g/L	10 g
<b>cQrex™ AAQ</b> N-Acetyl-L-alanyl-L-glutamine	Increases cell-specific productivity	>500g/L	5 g
cQrex™ AKG α-ketoglutaric acid	Titer-enhancing additive, protects cells against oxidative damage	>1000g/L	5 g
c <b>Qrex™ Carnosine</b> β-alanyl-L-histidine	Protects primary and stem cells against oxidative damage, influcences cell ageing	-	-

This information and all further technical advice are based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.

cQrex<sup>™</sup> – trademark of Evonik Industries AG and its subsidiaries.



### EVONIK NUTRITION & CARE GmbH

Health Care Business Line Pharma & Food Ingredients

cQrex@evonik.com www.evonik.com/cQrex