

Biosimilar Characterization for Regulatory Submission and Process Development

Overview of Company and Capabilities



Advancing Technologies and Applications of Proteome Analysis

Premises, Problem Statement and Solution

mAb / Vaccine / biotherapeutics development program
underway but

- Some or all components of physico-chemical characterization missing
- Product characterization needed for the purpose of regulatory filing and process development



Consider Shantani's more than a decade old
Expertise in Protein Analysis



Why Shantani

- International Scientific Team
- Cutting-Edge Technologies
- High-end Instrumentation

} Global Quality Standards

- Located in Pune, India
- Consultative Approach
- 'Fit-for-Purpose' Analysis


} Faster Turn Around Time

- Completed > 20 Biologics Characterization
- Supported > 12 Regulatory Submissions
- Worked with > 6 Biosimilar Developers

} Successful Track Record of Delivery



Drug Substance and Drug Product Characterization Workflow



- Intact Mass
- Peptide Mapping
- Glycan Profiling
- Charge Variant Analysis
- Circular Dichroism
- Fluorescence Spectroscopy
- Disulfide Bond Analysis
- FTIR
- AFM
- SEM / TEM
- Cytotoxicity Assays
- Host-Cell Protein / Host-Cell DNA / Aggregate Analysis



Biosimilar Characterization and Comparability Analysis



Intact Mass Analysis

- Typically carried out using ESI-Q-TOF (Agilent 6540 UHD)
- Standard and well-established protocols for Biotherapeutics
 - Underivatized / Derivatized
 - Glycosylated / De-glycosylated
 - Reduced / Non-reduced
- Successful looking at protein up to 240 kDa size





Determination of Intact Mass of mAbs using ESI-MS

Innovator's
Samples



Biosimilars



Primary mass matches, however, difference in glycosylated species observed !!!




Peptide-Map Analysis

- Standard and well-established protocols for single and multi protease, in-gel or in-solution, digestions of proteins
- Peptides are separated using RPLC and infused into ESI-QTOF set up for MS and MS/MS analysis
- Obtained masses are compared with the theoretical digest of the proteins using bioinformatic tools and peptide-maps are prepared
- Successfully prepared peptide maps of several (> 30 products) biotherapeutic proteins



Peptide Mapping of Biosimilar mAbs



Peptide Mapping of a Vaccine Product

- Mass data were compared against the theoretical digest of protein Uniprot ID P69899.
- Peptides containing Amino Acids in Red were identified.

MWRPSDSTVYVPPNPVSKVVATDAYVTRTNIFYHASSRLLAVGHPYFSIKRANKTVVPKVSGYQYRVFKVVLDPD
 NKFALPDSSLFDPTTQRLVWACTGLEVGRGQPLGVGVSGHPFLNKYDDVENSGSGGNPGQDNRVNVGMDYKQT
 QLCMVGCAPPLGEHWGKKGKQCTNTPVQAGDCPPLELITSVIQDGDMDVDTGFGAMNFADLQTNKSDVPIDICGTTT
 KYPDYLQMAADPYGDRLFFFLRKEQMFARHFFNRAGEVGEVPEVDTLIIKSGGNRTSVGSSIYVNTPSGSLVSSEAQ
 LFNKPYWLQKAQGHNNGICWGNQLFVTVDTRSTNMTLCASVTTSSSTYTNSDYKEYMRHVEEYDLQFIFQLCSIT
 LSAEVMAYIHTMNPSVLEDWNFGLSPPPNGTLEDYRYVQSQAITCQKPTPEKEKPDYKNLSFWEVNLKEKFSSE
 LDQYPLGRKFLQSGYRGRSSIRTGVKRPVASKASAAPKRKRAKTKR

(80% Coverage, 398 / 500 amino acids)



Secondary and Tertiary Structure Analysis

- If needed, samples are de-formulate using Gel-Filtration-Chromatography and then analyzed using Fluorescence and CD Spectroscopy
- Tryptophan Fluorescence Quenching experiments can be performed for in-depth analysis
- Successfully analyzed secondary and tertiary structures of several protein biotherapeutics of innovators and biosimilars





Steady State Fluorescence Spectroscopy

Innovator's Samples





Comparative Steady-State Fluorescence
Sample 5 to 7



Biosimilars




Comparative Steady-State Fluorescence
Sample 1 to 4



No Major
Deviation
Observed !!!

Circular Dichroism Studies



Sample #	Alpha-Helix	Beta-Sheets	Irregular
# 1	0%	53%	55%
# 2	0%	53%	55%
# 3	0%	53%	55%
# 4	0%	53%	55%
# 5	0%	53%	55%
# 6	0%	53%	55%
# 7	0%	53%	55%


Minor Variations in CD spectra observed, however, No significant Difference in secondary structure components !!!

Charge Variant Analysis

- Charge Variants are analyzed using HPLC-strong-cation exchange chromatography (SCX) and UV/Fluorescence Detectors
- Linear pH gradient is the key to success and well-optimized protocol is utilized for the purpose





HPLC-SCX-UV Probing of Charge Variants of mAb



Major difference in percentage contribution of Charge Variants Observed !!!

Disulfide Characterization in mAbs

Typically mAbs are digested with proteases in reduced and non-reduced conditions and obtained peptides are analyzed using LC-UV or LC-MS based workflows.




Glycan Profiling

- Standard protocols for N and O Glycan profiling
- Extracted N and O glycans are labelled with 2-AB and analyzed using LC-Flu based workflow and Glycans are identified using mass-spectrometry
- Fluorescence labelled glycans were also analyzed using Fluorophore Assisted Carbohydrate Electrophoresis (FACS)




N-Glycan Profiling of mAbs



Clear Variation in Mannose-5 composition in Biosimilar Samples were observed.

O-Glycan Profiling of mAbs



Species	Sample-1	Sample-2
Peeling	35%	36%
Gal	-	-
Neutral	15%	18%
Mono SA	43%	38%
Di SA	7%	7%

Species were confirmed by ABS and BTG digestion and Mass-Spectrometry.

No Major Difference Observed !!!


Lipid Analysis

- Mass-spectrometry and Thin-Layer Chromatography are utilized in characterizing the lipid components of the biotherapeutics
- Several lipid analysis method using mass-spectrometry has been standardized in our lab




Comparative Profile of Lipid Impurity in Low Molecular Weight Heparin Products (LMWH)

Lipid Mass Profile of Crude Heparin and Enoxaparin
Mirror Plot




Lipid Mass Profile of Lovenox and Enoxaparin
Mirror Plot



> 70 % Profiled lipid impurity were same between Enoxaparin (Biosimilar) and Lovenox (Innovator). By several supportive experiments it was concluded that lipid impurity profile of biosimilar and innovators product was same.

Our Focus = Protein Analysis



Our Key Strengths

- Globally Competitive Science
- High Ethical and Professional Standards
- Networked Operational Model for Cost-Effectiveness



Our Customer's Profile

East Coast-USA based one of the Largest Pharmaceutical Company

World's Largest Chemical Producer HQ - Germany



Glimpses of Shantani's R&D Center @ Innovation Park



Thank you

Connect for further discussions

Saba Naaz Siddique
saba@shantani.com

<http://www.shantani.com>



Advancing Technologies and Applications of Proteome Analysis

