

# Target Identification using Chemical Proteomics at Shantani

---

Chaitanya Saxena, Ph.D.  
CEO, Shantani

January, 2017



# Company Overview – Conception & Growth

## Chemical-Proteomics Based Biotechnology Company



# Premises, Problem Statement and Value

Discovery Program Identified a *bioactive* (small-molecule, peptide, protein, antibody) however the target(s), the cellular binding partners, of *bioactive* is not known

- Information of Target of *bioactive* will allow
  - Rational optimization of *bioactive*
  - Drug-Efficacy biomarker discovery, patient stratification and commercial differentiation by clarifying the action mechanism of *bioactive*
  - Early 'de-risking' of program by characterizing 'off-target' of the *bioactive*



# Shantani Identify Targets of Bioactive Molecules

Program's Target ID Need



Shantani Deploy (Appropriate Technology  
+ A Decade of Target Identification Expertise + Program  
Centric Business Model)



Deconvoluted Target Information




# Portfolio of Target-ID Technologies

Target Identification/Deconvolution is not Trivial =  
A single Tool / Technology May Not necessarily solve the problem for all

## Our Focus = Target ID

Our Expertise in the field allows us to evaluate the 'fit-for-purpose' technology and then we deploy appropriate Technology for right target



### Different Technologies

1) **UPT** = Unique Polymer Technology, 2) **SCLS** = Subcellular Location Specific Target Capture Technology, 3) **COMP** = 'in-silico' target ID workflows, 4) **TBB** = Traditional Bead/Biotinylated Molecule Based Method, 5) **TPP** = Thermal Proteome Profiling



# Comparative Profile of Target ID Technologies

Technolog(ies)	Key Technology Differentiator	Typical Deconvoluted Targets	Typical Time-Line (Weeks)	Comparative Cost
<b>Proprietary</b> Unique Polymer Technology (UPT)	'bait-molecule' derivatization not required, target enrichment based identification	8-10 (false positive rate ~40%)	2-3	\$\$\$
<b>Proprietary</b> Subcellular Location Specific Target Capture Technology (SCLS)	Target Capture from Live-cell in sub-cellular location specific manner	4-6 (false positive rate ~20%)	4-8	\$\$\$\$\$
<b>Proprietary</b> 'in-silico' Tools (COMP)	Robust, Fast and cost-effective	30-40 (false Positive rate ~80%)	1-2	\$
<b>Non-Proprietary</b> Bead/ Biotin Based Traditional Chemical Proteomics Technology (TBB)	Target Capture from cell-lysates	12-15 (false positive rate ~40%)	2-3	\$\$\$
<b>Non-Proprietary</b> Thermal-Proteome Profiling (TPP)	'bait-molecule' derivatization not required, target deconvolution base on bait-molecule induced thermal stabilization of target protein	10-12 (false positive rate ~50%)	3-5	\$\$\$\$



# Shantani's Proprietary Target-ID Technology Platforms

Technologies	Final Value
<p>Unique Polymer Technology (UPT) – Label Free Technology</p>	<ul style="list-style-type: none"> <li>• Quick Target Profiling allowing a 'Go/No Go' development decision for phenotypically screened compounds</li> <li>• Rapid target profiling of multiple compounds to save time and cost</li> </ul>
<p>SubCellularLocation Specific Target Capture Technology (SCLS) – Labelled Technology</p>	<ul style="list-style-type: none"> <li>• Precise Target Information driving rational lead development</li> <li>• Low False Positives saving time during target validation</li> </ul>
<p>Computational Supreme (Comp-S) – Label Free Technology</p>	<ul style="list-style-type: none"> <li>• Significantly Narrowed down list of targets to be validated, saving time</li> <li>• Assist in bringing selectivity while lead development / optimization</li> </ul>



# Unique Polymer Technology(UPT)

---





# Unique Polymer Technology (UPT)

## Key Advantages

- Derivatization of test-molecule is not needed
- Target Deconvolution can be completed within 2-3 weeks
- Low False Positive Identification Rate = Faster Validation




# Immobilization without derivatization : Working Hypothesis

- Weak molecular interaction forces of a organic molecule can be used in immobilizing them on a surface that provides complementary weak interactions
- ↓
- Sum of multiple weak interactions is strong enough to allow the molecule to stay (immobilized) on the surface
- ↓
- The molecule can stay on the surface for long enough they can be used as molecule specific matrix for affinity capture of protein target



# Design of Unique Polymer for Chemical Proteomics




**Multiple Possible non-covalent Affinity Binding Interactions for small-molecules**

Ionic, Hydrophobic / Hydrophilic  
Vander-waals, pi-pi, cation-pi, H-bonds

**Optimized Surface Angel Properties for attachment of polymer to**

Glass and/or Plastic in different dissolution Phase

**Overall Amphiphilic in Nature**



 = Unique Polymer



# Preparation of Test-Molecule Specific Affinity Matrix

Test molecules are dissolved in appropriate solvent



Test-molecules are layered on the polymeric matrix and allowed to stay for appropriate time for immobilization



Matrix is washed extensively and amount of test-molecule in washes confirms the extent of immobilization




Example: Polymer layered in a Eppendorf tube and later a bisindolylmaleimide compound was immobilized



# Retention Behavior of a few Small-Molecules on Unique Polymer

Retention Behaviour of the Compounds on Polymeric Surface  
During Washing of Compound Coated Polymers



# Molecule Orientation on Polymer is critical

- Which interaction force is dominant and whether molecules are oriented in linear or random fashion ?
- If orientation is not random then a particular site of molecule that is responsible for its interaction with the target protein may not be exposed !!!

Let's consider one major 'weak-interaction' force at a time




# Orientation of Molecule is Random in Nature

## Let's consider Only Hydrophobic Interactions

If only hydrophobic interactions are playing a role in compound immobilization then different compound of about similar logP values should be retained similarly

Retention Behaviour of the Compounds on Polymeric Surface During Washing of Compound Coated Polymers



Retention behavior of compounds having about similar logP value is significantly different

= **More than one type of interaction forces are working together** to retain the molecule on the polymer


= Random orientation of the molecule on polymer is plausible factor for multiple weak interactions



# Orientation of Molecule is Random in Nature

## Let's consider Only Ionic Interactions

If only ionic interactions are playing a role in compound immobilization then polymer coating and washing with buffer having different pH will ionize the compound differently and retention will be affected



Retention behavior of same compound coated and washed at different pH did not significantly change its retention behavior

= More than one type of interaction forces are working together to retain the molecule on the polymer

= Random orientation of the molecule on polymer is plausible factor for multiple weak interactions






# How much amount of Molecule is needed on the polymer for good capture

- Molecules are retained but how much amount is needed on the polymer ?
- Aqueous buffers are used to characterize the retention, what if the incubation with cell-lysate takes the molecule '*off-the-polymer*' ?


Let's Answer it by looking at Qualification Criteria of Matrix Preparation and Target Capture 



# Method Qualification *is a function of* Amount of Molecule Immobilized on given surface area of Polymer)



# Unique Polymer Technology (UPT) Work-Flow



Reference: *Shantani's Proprietary Technology.*  
 Patent Application :PCT/IN2017/000002



Capture of GSK3 protein, a well-established  
protein target ( $K_d \sim 19$  nM) of  
Bisindolylmaleimide-III using UPT



# Capture of Known Target of Bis-III using Polymer Layered in 1.5 ml tube – *Very First Capture Experiment using UPT*

**Step-1)** Affinity Matrix of underivatized molecule




Polymeric Affinity Matrix of Bis-III

Control Polymeric Matrix



**Step-2)** Affinity Chromatography and SDS-PAGE and Western-blot analysis of Known Target Protein




Primary Target was Specifically Captured but capture efficiency was not great !!!




# Capture of Target of Bis-III using Polymer Layered on Slide/CHIPs

**Step-1)** Affinity Matrix of underivatized molecule



**Step-2)** Affinity Chromatography and SDS-PAGE and Western-blot analysis of Known Target Protein



Primary Target was Specifically Captured and Amount of captured protein was significantly higher.



# Other Identified Targets of Bis-III using Polymer Layered on Slide/CHIPs

Uniprot_ID	Protein Description	Maximum Number of Unique Peptides Identified	Q-Value (%)
Q13418	Integrin-linked protein kinase	7	0
Q70UQ0	Inhibitor of nuclear factor kappa-B kinase	6	0
P28482	Mitogen-activated protein kinase 1	6	0
P60891	Ribose-phosphate pyrophosphokinase 1	5	0
E9PF82	Calcium/calmodulin-dependent protein kinase type II	4	0
P49841-2	Glycogen synthase kinase-3 beta	3	0
P63208	S-phase kinase-associated protein	3	0
P51570-2	Galactokinase	3	0



# Capture of protein targets of new compounds – Case Studies






# Case Study # 1: Capture of Target(s) of Phenotypic Hit (Research Institute, Seattle)

- **Scenario:** Antibacterial compound not working through bacterial cell-wall disruption was identified. Target not known. Molecule can not be developed further.
- **Client Provided: Molecule** - 15 mg | **Cell-Lines** – 9 mg Bacterial Lysate
- **Time-lines at Shantani:** 1 month
- **Final Deliverable:** Very small list of highly putative targets of molecule.

**Step-1)** Preparation of Affinity Matrix of un-derivatized molecule & Control Matrix



**Step-2)** Affinity Chromatography and SDS-PAGE Analysis




Comparative profile of proteins obtained from control and drug loaded matrix on SDS-PAGE does not appear distinct.



Mass-Spec based identification of all proteins will tell the right story.



# Case Study # 1: Target Deconvolution Process




# Case Study # 1: Final Deconvoluted Target

Uniprot_ID	Protein Description	Maximum Number of Unique Peptides Identified	Protein Sequence Coverage (%)	Q-Value (%)
P9WQP1	Diacylglycerol acyltransferase/mycolyltransferase Ag85B	2	8.92	0
P9WN21	Fructose-1,6-bisphosphatase class 2	3	13.81	0
P9WFBV1	Leucine-tRNA ligase	3	3.82	0
I6Y7V6	Acyl-CoA ligase FadD31	2	4.52	0
P9WMJ7	Chaperone protein HtpG	2	5.87	0
I6YBZ8	3-hydroxyacyl-thioester dehydratase	6	28.97	0
P9WHV1	Gamma-glutamyl phosphate reductase	8	27.71	0
P9WMJ9	Chaperone protein DnaK	31	64.00	0




# Case Study # 2: Capture of Target(s) of Phenotypic Hit (Research Institute, EU)

- **Scenario:** Molecule disrupts Notch signaling pathway, hard to derivatize, upon derivatization loosed activity, target not know.
- **Client Provided: Molecule - 12 mg | Cell-Lines - RPMI-8402 & Cell Lysate: 9 mg**
- **Time-lines at Shantani: 2 Weeks**
- **Final Deliverable:** Very small list of highly putative targets of molecule.



# Case Study # 2: Results

**Step-3) SDS-PAGE Analysis**



→ **Step-4) Target ID using Mass-Spectrometry**

→ **Step-5) Target-Deconvolution**

Total Proteins Identified = 253

↓  
Deconvoluted Specific Targets = 6

Protein Class / Name	Specificity Ratio
Methyl Transferase	1
Fibronectin Binding Protein	1
Mitochondrial Transport Protein	1
GTPase activating Protein	1
Adapter Protein for T-Cell Signalling	1

A novel non-HDAC epigenetic target was outlined and later validated by client.



# Case Study # 3: Target Capture to explain the phenotype (Client – a large Indian Corporate)


- **Scenario:**
  - Molecule and its targets are established in non-human systems. In human cell-lines molecule significantly inhibits PGD2 production and increases proliferation, a property of interest to client.
  - Protein target(s) of the molecule in human system are not known and hence program can not be developed further in rational manner.
- **Client Provided:**
  - **Molecule** - 50 mg (used 12 mg) | **Cell-Lines** – CONFIDENTIAL
- **Time-lines at Shantani:** 5 Weeks
- **Final Deliverable:** Very small list of highly putative targets of molecule.




## Case Study # 3:

### Step-1) Preparation and Qualification of Test Molecule Specific Affinity Matrix

Project # PS0035, Compound # CONFIDENTIAL  
Retention Pattern Matrix-1



Project # PS0035, Compound # CONFIDENTIAL  
Retention Pattern Matrix-2



Higher amount of molecule was captured on Matrix-1 and was chosen for further experimentation.



## Case Study # 3:

**Step-2, 3 & 4) Affinity Chromatography, Target Deconvolution and Prioritization of Target Validation Experiments**

5 Target were identified.



Computational docking studies was utilized to further prioritize the identified putative targets for target validation efforts.



List was narrowed down to 3 targets

Protein Class / Name	PDB_ID	Predicted Ki ( $\mu\text{M}$ )	Relative Docking Rank & Confidence
Synaptobrevin homolog	3KYQ	3.93	90%
ARL1	1HUR	30.90	30%
TRPV2	2ETA	86.49	33%
Exportin-T	3IBV	81.42	18%
huRNP	2OT8	18.53	15%





## Case Study 3 #:


One of the Prioritized Target is ARL1 (ADP-ribosylation factor like protein 1)



ARL1 shares 57% of amino-acid sequence with ARF1 (ADP-ribosylation factor 1)



ARL1 is directly implicated in PGD2 pathway



Identified Targets allowed further development of the program.




# Case Study # 4 Toxicity Profiling (Client – one of the biggest global biopharma company)

- **Scenario:**
  - BACE-I inhibitor, one of the promising therapies for Alzheimer's
  - Preclinical compounds shows ocular toxicity in animal models
- **Client Provided:**
  - **Molecules** - 20 mg (Multiple molecules) | **Cell-Lines** – RPE (relevant to Ocular Toxicity)
- **Time-lines at Shantani:** 5 Weeks
- **Final Deliverable:** Very small list of highly putative targets of molecule.



# Case Study - 4

## Step-1) Preparation and Qualification of Test Molecule Specific Affinity Matrix



## Case Study - 4

### Step-2, 3 & 4) Affinity Chromatography and Target Deconvolution

12 Targets but belonging to only 2 different class of protein were identified



Relevance of both the protein classes with Occular Toxicity was evaluated



List was narrowed down to 3 targets



One of the Cathepsin Family Member was identified as primarily responsible for the toxicity

**Value Added** – Could pinpoint target protein responsible for the toxicity, 'lead-compound' selectivity process got a significant boost.




# Conclusions from UPT based Chemical-Proteomics Methodologies

- 'Bait-molecule' derivation and SAR information not needed for Target ID
- Target Deconvolution can be carried out as fast as in 2 weeks
- False-positive identification rates though appears higher can be controlled by running multiple experimental replicates
- UPT can be effectively used for large and fast screening of small-molecule targets.




# Several Examples


Retention Compound




Re




Re



Retention P Compound




Compound # LS 5968094 / Matrix # 4




Retention Pattern of H2-Gamendazole


Retention Pat on Diff




Retention on Di




on



Retention Pat




% compound remain retained on the polymer matrix after different washes



Matrix # 1  
Matrix # 2



# Target ID and Client Engagement Workflow



# SubCellular Location Specific Target-Capture Methodologies (SCLS)

---





# Hypothesis

## Target-Capturing Probes that Can

- Identify sub-cellular compartment of molecular activity
- Be used at functionally relevant concentration for target capture
- Be Recovered from live cells

Will allow capture of rightful targets of small-molecule from physiologically relevant live-cell systems



# Step-1: Identification of Sub-Cellular Location


a) Location specific peptide coupled "bait molecule"

b) Biological System

c) Activity Read-Outs




d) Subproteome Selection

Membrane Probe




Membrane

Cytoplasmic Probe



Cytoplasm


Nuclear Probe




Nucleus



# Step-2: Target Capture from Sub-Cellular Location and Target Identification




# Background: Development of Sub-cellular Location Specific Peptides



# Properties of Peptides

- Penetrate wide variety of mammalian cells and remain confined in sub-cellular location specific manner
- Not toxic to the cells
- No known peptidase sequence - Do not degrade in the cell
- 6-10 amino-acids long - Do not form secondary structures






# SCLS –Validation Experiments

---



# Proof of Concept Experiments (PoCs)


Three known molecules with known targets and sub-cellular target location were chosen for PoC Experiments

			
<b>Molecule</b>	<b>Bis-III</b>	<b>BMS345541</b>	<b>Histamine</b>
Target	GSK3-beta	IKK-2	H1-Receptor
Location	Nucleus / Cytoplasm	Cytoplasm	Membrane
Functional Activity	Cytotoxic	Inhibition of LPS Induced Nitric –Oxide Production	Increase Intracellular Calcium




# Preparation of Bis-III Probes

**Step-1:** Analysis of Structure Activity Relationship (SAR) for identifying site for peptide coupling



**Step-2:** Coupling of Location-Specific Peptide with Bis-III followed by HPLC based purification and Mass-Spectrometry Based Characterization of Probe






# Sub-cellular Location Specific Functional Activity of Bis-III Probes

Bis-III was coupled to three sub-cellular location specific proprietary peptides



Cytotoxic Assay in HeLa Cells



Bis-compound and all probes and control probes at 80  $\mu$ M concentration.


Bis-III probes that were targeted to the nucleus exhibited similar activity as the free Bis-III



Target located in nucleus is responsible for the activity




# Capture of Target from Nuclear Fractions



GSK3-beta ⇒

Probes = Coupled with Bis-III




Target was specifically captured with Probe.



# Sub-cellular Location Specific Functional Activity of BMS345541 Probes

BMS345541 was coupled to three sub-cellular location specific proprietary peptides

**Functional Assay:** Inhibition of LPS induced Nitric Oxide Production in RAW cells




- Basal Control = No induction of Nitric Oxide by LPS
- BMS compound = 12.5 µM
- All peptide probes and control probes = 25 µM

BMS probes that were targeted to the cytoplasm exhibited similar activity as the free BMS compound


Target located in cytoplasm is responsible for the activity



# Capture of Target from Cytoplasmic Fractions



Probes = Coupled with BMS345541




Target was specifically enriched with Probe.



# Sub-cellular Location Specific Functional Activity of Histamine Probes

Histamine was coupled to three sub-cellular location specific proprietary peptides

**Functional Assay:** Increase in intracellular calcium in HeLa Cells






Histamine-compound and all probes at 200  $\mu$ M concentration.

All the histamine probes showed the activity

Every peptide probe shows the activity, i.e. target is located on the membrane because every probe will either hit or cross the membrane activating increase in calcium



# Capture of Target from Cytoplasmic Fractions



Target was specifically captured with Probe.



**SCLS**  
**VS**  
**Traditional Bead Based Method**

---



# Direct Comparison: SCLS vs Traditional Bead Based Method – Molecule # 1

Bait-molecule = Bis-III

## Workflow Comparison

	SCLS	Bead Based
<b>Biological System</b>	Physiologically Relevant Live-cells	Cell-Lysate
<b>Target Capturing Probe</b>	Nuclear Specific Probe of Bis-III (based on the functional activity)	Bis-III immobilized on sepharose 6B Epoxy beads
<b>Affinity Interaction</b>	55 $\mu$ M Bis-III Probe in live cell	25 $\mu$ l Bis-III on beads with cell lysate
<b>Affinity Pull-Down</b>	Antibody Against the Probe	Bis-III on beads
<b>Protein Analysis</b>	SDS-PAGE followed by in-gel trypsin digestion and Mass-spec based protein identification	SDS-PAGE followed by in-gel trypsin digestion and Mass-spec based protein identification

Results on next slide





# Direct Comparison: SCLS vs Traditional Bead Based Method – Molecule # 1

Criteria	SCLS	Traditional Bead Based
Primary Target Identified	YES	YES
Confidence in Target Identity	Strong (4 unique peptides identified)	Weak (only 1 unique peptide identified)
Total Number of Deconvoluted Protein Targets	15	16
Capture from Native Cellular Environment	YES	NO
Capture using Functionally Relevant Ligand	YES	NO
Number of Target(s) (with Known Strong Affinity to Molecule) Identified	4	1
Number of Target(s) (with Known Low Affinity to Molecule) Identified	2	3
Dynamics Range of Capture	Medium (Primarily identify high affinity binding partner)	High (Identify high & low affinity binding partners)
False Positive Capture	Low	High

SCLS provides better target deconvolution.



# Direct Comparison: SCLS vs Traditional Bead Based Method – Molecule # 2

Bait-molecule = SB202190, Known Target = p38

## Workflow Comparison

	SCLS	Bead Based
<b>Biological System</b>	Physiologically Relevant Live-cells	Cell-Lysate
<b>Target Capturing Probe</b>	Cytoplasm Specific Probe of SB202190 (based on the functional activity)	SB202190 immobilized on sepharose 6B Epoxy beads
<b>Affinity Interaction</b>	25 $\mu$ M Bis-III Probe in live cell	25 $\mu$ l Bis-III on beads with cell lysate
<b>Affinity Pull-Down</b>	Antibody Against the Probe	Bis-III on beads
<b>Protein Analysis</b>	SDS-PAGE followed by in-gel trypsin digestion and Mass-spec based protein identification	SDS-PAGE followed by in-gel trypsin digestion and Mass-spec based protein identification

Results on next slide



# Targets Identified for SB202190

Protein Annotation	No. of Unique Peptides Identified	Sequest Xc For highest Matched peptide	Tandem e value For highest Matched peptide	Specificity Ratio
Mitogen-activated protein kinase 14 isoform 2 (p38)	9	5.51	0.000	1.00
Isoform 1 of Glycogen synthase kinase-3 beta (GSK3-beta)	5	5.49	0.000	1.00
Lactoylglutathione lyase (GLO1)	5	5.31	0.003	1.00
Isoform 1 of Casein Kinase I delta (CKId)	3	4.35	0.000	1.00
STE20/SPS1-related proline-alanine-rich protein kinase (STK39)	2	6.13	0.000	1.00
Isoform Alpha-2 of Mitogen-activated protein kinase 9 (JNK2)	5	6.54	0.005	0.83

Robust Target Identification



# Direct Comparison: SCLS vs Traditional Bead Based Method – Molecule # 2

	SCLS	Bead Based Method
Number of Deconvoluted Target	6	18
Validation Time (assuming one target / month, validation at in-vitro/cellular level only)	6 months	18 months

Better Target Deconvolution = Faster Validation =  
Program Go/No-Go Decision = Time and Cost Saving



# Validation of Identified Target

Protein Annotation	Interaction Efficiency
Mitogen-activated protein kinase 14 isoform 2 (p38)	Inhibition 95.95% (Kinase Panel) IC50 = 38 nm
Isoform 1 of Glycogen synthase kinase-3 beta (GSK3-beta)	Inhibition 68.65% (Kinase Panel)
Lactoylglutathione lyase (GLO1)	ND
Isoform 1 of Casein Kinase I delta (CKId)	Inhibition 88.72% (Kinase Panel)
STE20/SPS1-related proline-alanine-rich protein kinase (STK39)	Inhibition 63.25% (Kinase Panel)
Isoform Alpha-2 of Mitogen-activated protein kinase 9 (JNK2)	Inhibition 72.25% (Kinase Panel)

Identified Targets were True Positive Targets of SB202190




# From Quantitative Target Information to Action Mechanism and Tox Prediction

For SB202190 compound the identified proteins were rank-ordered based on their abundance and interaction ability and were mapped on to the canonical pathways to predict that molecule may induce apoptosis

## Pathway Knowledge Base

Identified Proteins if inhibited together will lead to Apoptosis.



# SB202190 Does Induce Apoptosis

## A selective inhibitor of p38 MAP kinase, SB202190, induced ...

[www.sciencedirect.com/science/article/pii/S0925443900000454](http://www.sciencedirect.com/science/article/pii/S0925443900000454) ▼

by H Karahashi - 2000 - Cited by 36 - Related articles

A selective p38 MAP kinase (p38 MAPK) inhibitor, **SB202190**, induced apoptotic cell death of a macrophage-like cell line, J774.1, in the presence of ...

## Induction of Apoptosis by SB202190 through Inhibition of ...

[www.jbc.org/content/273/26/16415.full.pdf](http://www.jbc.org/content/273/26/16415.full.pdf) ▼

by S Nemoto - 1998 - Cited by 279 - Related articles

... addition, SB202190 was able to potentiate apoptosis induced by Fas(APO-1) ...  
These re- sults indicate that **SB202190 induces apoptosis** through activation of ...

## Induction of Apoptosis by SB202190 through Inhibition of ...

[www.jbc.org/content/273/26/16415.full](http://www.jbc.org/content/273/26/16415.full) ▼

by S Nemoto - 1998 - Cited by 279 - Related articles

Jun 26, 1998 - In contrast, expression of p38 $\alpha$  induced cell death mildly. These results indicate that **SB202190 induces apoptosis** through activation of ...

## A selective inhibitor of p38 MAP kinase, SB202190, induced ...

[www.ncbi.nlm.nih.gov/pubmed/11040446](http://www.ncbi.nlm.nih.gov/pubmed/11040446) ▼

by H Karahashi - 2000 - Cited by 36 - Related articles

Oct 18, 2000 - A selective inhibitor of p38 MAP kinase, **SB202190**, induced apoptotic cell death of a lipopolysaccharide-treated macrophage-like cell line, ...

## Induction of apoptosis by SB202190 through inhibition of ...

[www.ncbi.nlm.nih.gov/pubmed/9632706](http://www.ncbi.nlm.nih.gov/pubmed/9632706) ▼

by S Nemoto - 1998 - Cited by 279 - Related articles

J Biol Chem. 1998 Jun 26;273(26):16415-20. **Induction of apoptosis by SB202190** through inhibition of p38 $\beta$  mitogen-activated protein kinase. Nemoto S(1) ...



# Computation Workflow

---






# Computational Supreme (Comp-S)

## Key Advantages

- Utilizes Structural and Functional Features of 'test-molecule' when comparing historical database
- Quick and cost-effective
- Low False Positive Identification Rate = Faster Validation




# Case # 5: Deconvolution of Targets of a Few Known Small Molecules using Shantani's COMP workflow



Bisindolylmaleimide-III

	# of Target Identified	# of Actual Target (Kd < 500 nm)	Success Ratio (Total / Actual Target)
Algorithm-1	26	1	4%
Algorithm-2	4	2	50%
Algorithm-3	16	2	13%
<b>Shantani-Algorithm</b>	<b>2</b>	<b>1</b> (Identified Target =PKC- $\alpha$ )	<b>50%</b>



BMS-345541

	# of Target Identified	# of Actual Target (Kd < 500 nm)	Success Ratio (Total / Actual Target)
Algorithm-1	17	1	6%
Algorithm-2	8	2	25%
Algorithm-3	19	2	11%
<b>Shantani-Algorithm</b>	<b>2</b>	<b>1</b> (Identified Target =IKK- $\beta$ )	<b>50%</b>

## Eventual Value and 'GO' Decision

Shantani's COMP workflow allows identification of a few but rightful targets



# Identification of Target of Novel Molecule

---



# Identification and Validation of Target and Understanding MoA of a Novel Imidazoline

**Molecule:** BL11282 Stimulates insulin secretion in unknown fashion.

- a) BL11282 does not block ATP-dependent K<sup>+</sup> channels
  - b) BL11282 activity appears to be sensitive to the inhibition of PKA
  - c) The increase in intracellular calcium upon BL11282 administration makes only a minor contribution to insulin secretion
  - d) BL11282 directly influences exocytosis process.....
- The Cellular Binding Partners of BL11282 are not known and its action mechanism is elusive hence the molecule can not be rationally optimized...




**BL11282**

**In this PoC we will utilize our chemical-proteomics tools and strategies to identify and validate primary targets of BL11282**



# Synthesis of BL11282 derivative for preparation of target capturing probe




References: BL11282 - Structure-Activity Relationship (SAR) data

- Bioorganic and Medicinal Chemistry 15: 3284-3265 (2007)
- Bioorganic and Medicinal Chemistry 15: 6782-6795 (2007)



# Functional Activity of BL11282 derivatives and Target Capturing Probes



Peptide 1, 2 and 3 are three different peptides that targets three different locations of the cell.




BL11282 was coupled to these peptides and functional activity (insulin secretion ability) of the constructs was evaluated.


Sub-Cellular Location of the Target was Identified



# Utilization of Traditional Bead Based Method for Target Pull-Down from Specific Cellular Fraction




# Distinct Elution Profile over Wash-Condition Was Optimized





# Distinct Protein Profile Over Control Experiments

CBB Stained Gel-Image



Lane 0 = Std

Lane 2 = Control Matrix

Lane 4 = Elution from BL11282 Matrix

Lane 6 = Elution from BL11282 Matrix  
(cell-lysate pre-incubated with BL11282)

Lane 1,3, 5, 7 & 8 = Buffer

\*Each of these lanes are proteins pooled from three different experiments and at least two such runs were carried out for protein identification. Protein bands were cut, in-gel trypsin digested and identified using mass-spectrometry based workflow.

**Note:**The proteins were separated only for 7 minutes on SDS-PAGE gels.




# Deconvoluted Proteins were Subjected to 'in-silico' docking analysis with BL11282

Small-molecule affinity based chromatographic methods will elute **protein complexes associated** with the target protein rather than **The Only Target Protein.**



Proteins that may be directly interacting with the small-molecule can be prioritized based on their binding probability to the targets




# Prioritization of Proteins using Predicted Free-Energy of Interaction

Top 10 Target Proteins	Predicted Free Energy (kcal/mol)	Predicted (Kd) $\mu\text{M}$
1	-8.71	0.41
2	-8.61	0.49
3	-8.23	0.93
4	-8.03	1.30
5	-7.94	1.50
6	-7.54	2.96
7	-7.52	3.07
8	-7.43	3.60
9	-7.44	3.77
10	-7.4	3.79




# BL11282 Inhibits Target 1 and 2 – Biophysical Validation




Inhibition of Target – 1 by BL11282  
(Spectrophotometric Assay Monitoring Product of the  
Enzyme Reaction)



Inhibition of Target – 2 by BL11282  
(ELISA based Cell Based Assay)



# siRNA and Pharmacological Inhibition Based Functional Validation of Target(s)




- Cell Type – MIN6 (Mouse Insulinoma)
- Target knock down by siRNA was optimized and established separately.
- Data are plotted after subtracting respective control.

Both the targets if inhibited together increases insulin secretion !!!



# Conclusions From BL11282-Chemical Proteomics

- Study of molecule SAR allowed preparation of functionally active target capturing probe
- Appropriate chemical proteomics tools were used to identify cellular binding partners of BL11282
- Identified binding partners were thoroughly analyzed using in-silico tools and their putative relation with phenotype to deconvolute the targets
- Prioritized targets were validated using biophysical and molecular biology tools



# Shantani Identify Targets and MoA of Bioactive Molecules

**Program's Target ID Need**



**Shantani Deploy (Appropriate Technology + A Decade of Target Identification Expertise + Program Centric Business Model)**



**Deconvoluted Target Information**



# Our Key Strengths

- Globally Competitive Science
- High Ethical and Professional Standards
- Networked Operational Model for Cost-Effectiveness





# Our Customer's Profile

East Coast-USA based  
one of the Largest  
Pharmaceutical  
Company

World's Largest  
Chemical Producer  
HQ - Germany



THE UNIVERSITY OF TEXAS  
**MD Anderson**  
~~Cancer~~ Center



**KU** MEDICAL  
CENTER  
The University of Kansas

**reclaimrx**  
recover and regenerate

**EPFL**  
ÉCOLE POLYTECHNIQUE  
FÉDÉRALE DE LAUSANNE



**praj**  
Innovate • Integrate • Deliver

**national centre for cell science**  
An Autonomous Institution of the Department of Biotechnology, Govt. of India



**BAREFEET**  
Analytics

**SHIV NADAR UNIVERSITY**



**BIONEEDS INDIA PRIVATE LIMITED**  
INTEGRATED DISCOVERY, DEVELOPMENT AND REGULATORY SERVICES

**SEAGULL**  
BioSolutions Pvt. Ltd.

# Shantani R&D Center @ Innovation Park



# Thank You.

Connect for further discussions

Chaitanya Saxena, Ph.D.

csaxena@shantani.com

+91-9975447489, +91-20-64103918

<http://www.shantani.com>



Advancing Technologies and Applications of Proteome Analysis

