
Functional polymers to take control of your release profile

Versatility and reliability for oral solid dosage forms

Eudragit®





TAKE CONTROL OF YOUR RELEASE PROFILE WITH EUDRAGIT®

IMMEDIATE RELEASE:

- Taste and odor masking, light and moisture protection
- · Swallowability improvement
- Soluble in gastric fluid up to pH 5.0
- · Swellable and permeable above pH 5.0
- · Low viscosity
- · High pigment binding capacity
- · Excellent adhesion
- High effectiveness in thin coatings

EUDRAGIT® E 100 Granules

EUDRAGIT® E 12,5 Organic solution

EUDRAGIT® E PO Powder
EUDRAGIT® E PO ReadyMix Powder
EUDRAGIT® IR ReadyMix White Powder

DELAYED RELEASE:

- Easily combined to achieve a specific dissolution pH
- · Highly stable, with rapid dissolution for precise targeting
- Well-defined characteristics to safeguard the gastric mucosa from aggressive actives
- Able to protect the transit of APIs sensitive to gastric fluid

Small Intestine:

EUDRAGIT® L 30 D-55 Aqueous dispersion soluble above pH 5.5 EUDRAGIT® L 100-55 soluble above pH 5.5 Powder EUDRAGIT® FL 30 D-55 Aqueous dispersion soluble above pH 5.5 EUDRAGIT® L 100 Powder soluble above pH 6.0 EUDRAGIT® L 12,5 soluble above pH 6.0 Organic solution EUDRACAP™ enteric NEW Capsule soluble above pH 5.5

EUDRACAPTM Select NEW Capsule customizable

Value-adding products:

PlasACRYL® HTP20

Easy-to-use glidant and plasticizer premix, specifically designed for

EUDRAGIT® L 30 D-55 formulations

Colonic delivery:

EUDRAGIT® S 100	Powder	soluble above pH 7.0
EUDRAGIT® S 12,5	Organic solution	soluble above pH 7.0
EUDRAGIT® FS 30 D	Aqueous dispersion	soluble above pH 7.0
EUDRAGIT® FS 100	Powder	soluble above pH 7.0


Value-adding products


PlasACRYL® T20

Easy-to-use glidant and plasticizer premix, specifically designed for

EUDRAGIT® FS 30 D formulations

Monographs and DMFs are available for all EUDRAGIT® polymers upon request.





SUSTAINED RELEASE

- Multiple combination options to control delivery throughout the GI tract
- Solutions for once-per-day dosing, pulsatile release and other custom profiles
- Insoluble with pH-independent permeability
- · Flexible options for high or low permeability

EUDRAGIT® RL PO Powder EUDRAGIT® RL 100 Granules

EUDRAGIT® RL 30 D Aqueous dispersion EUDRAGIT® RL 12,5 Organic solution

- · Insoluble, high permeability, pH-independent swelling
- Customized release profiles by combining with EUDRAGIT® RS at different ratios

EUDRAGIT® RS PO Powder EUDRAGIT® RS 100 Granules

EUDRAGIT® RS 30 D Aqueous dispersion EUDRAGIT® RS 12,5 Organic solution

- Insoluble, low permeability, pH-independent swelling
- Customized release profiles by combining with EUDRAGIT® RL at different ratios

EUDRAGIT® NM 30 D Aqueous dispersion

- · Insoluble, low permeability, pH-independent swelling
- · Highly flexible, suitable as matrix formers

SOLUBILITY AND BIOAVAILABILITY ENHANCEMENT

- Polymer, process technology and formulation solutions
- Specialized in amorphous solid dispersions
- · Broad miscibility with APIs and excipients
- Expedited drug development through miniaturized screening tool
- Well-defined characteristics for spray drying, hot melt extrusion and top-spray granulation
- · Down-streaming and further processing options

EUDRAGIT® E POEUDRAGIT® L 100EUDRAGIT® E 100EUDRAGIT® S 100EUDRAGIT® L 100-55EUDRAGIT® FS 100

Chematic® for equipment cleaning

- Customized detergents for EUDRAGIT® RL, RS, NE, NM, L, FL, FS and S
- · Suitable for WIP, CIP and manual processes
- · Highest cleaning effectivity
- Fast and safe cleaning process

This information and all further technical advice are based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.



EVONIK OPERATIONS GMBH

Health Care Business Line Drug Delivery & Medical Device Solutions

eudragit@evonik.com www.evonik.com/eudragit