



CONTAINMENT & SAMPLING TECHNOLOGY

- Bespoke gloveboxes to allow for contained small scale processing of HAPI compounds
- Gloveboxes may be linked directly to existing or new fume cupboards via material transfer airlocks to segregate liquid processing within the fume cupboards and solids handling within the glovebox
- Gloveboxes may contain a variety of different processing equipment, such as tray dryers, Nutsche filters, rotary evaporators, small scale reactors and crystallisers, etc.
- Dispensing and packaging of final solids compound may be carried out within the glovebox



ProSys Laboratory Small Scale Isolators
ProSys Sampling Systems Ltd. • info@prosys.ie • www.prosys.ie • +353 21 4613890



www.prosys.ie