



CONTAINMENT & SAMPLING TECHNOLOGY

- Pharmaceutical active ingredients dispensing stations
- Can be designed to achieve a range of OEL levels
- Fully automated nitrogen purge system and pressure control circuit
- Designed to dispense either into continuous liner, docked drums or into Split Butterfly Valve (integrated into the base of the glovebox)
- Drum docking to the base of the glovebox for either the removal or addition of small quantities
- Integrated weighscales



ProSys Dispensing and Subdivision Isolators

• ✉ info@prosys.ie • 🌐 www.prosys.ie • 📞 +353 21 4613890



www.prosys.ie