




CONTAINMENT & SAMPLING TECHNOLOGY

• info@prosys.ie • www.prosys.ie • +353 (0)21 4853900

ProSys Recirculation Top Reactor Sampler



- Draws a representative sample from the top of a reactor vessel
- Handles slurries and crystals up to 5mm diameter
- PAT integration : pH measurement, conductivity etc.
- Can be integrated in a mobile skid for use on multiple reactors
- PTFE tubular diaphragm pump design allows 3m dry suction lift.
- Allows small quantities solids charge (seeding), can be integrated with a glovebox containment system.



www.prosys.ie