



CONTAINMENT & SAMPLING TECHNOLOGY

- Dry chamber glovebox for analytical QA facility
- Achieves less than 1% RH by continuous Dew Point analysis
- Fully integrated recirculation desiccator system
- Nitrogen purge facility for airlock chamber and working chamber
- Air purge system to remove all nitrogen from either chamber before windows are unlocked for opening
- Full touch screen HMI that allows for monitoring the running and maintenance status of the isolator as well as fault diagnostics
- 'Traffic Light' type indicator to allow quick status check of isolator
- Integrated power sockets, service hook up points and data transfer points



ProSys Dry Boxes

ProSys Sampling Systems Ltd. • info@prosys.ie • www.prosys.ie • +353 21 4613890



Cert No. 3342

www.prosys.ie